

Analytical Sample Preparation: Flexible, Easy and Rugged

Increase efficiency, productivity and meet highly demanding quality standards with our robotic platform for **faster and better analytical sample preparation**

The **SWING SP** is a modular robotic platform enabling automated and traceable substance dispensing (even while processing). Beyond the simple mixing of solids and liquids, the **SWING SP** performs more complex tasks like viscous material handling, filtration, labelling, barcode tracking and integrates with analytical instruments to provide a smooth and complete workflow.





SWING SP robotic platform

It has been designed with flexibility in mind, allowing you to reconfigure the platform and add or replace modules as your needs evolve.

Any type of chemical (liquids, solids, powders, granulates, spheres, viscous materials, waxes etc.) can be gravimetrically dispensed into the disposable vial/reactor format of your choice. Any sample preparation workflow can be managed by the robotic platform by adding and/or changing functionality as required (choice of up to 40 tool-features).


Workflow example:



Click on the picture



Click on the picture



Click on the picture



Click on the picture

Challenge us with your most difficult workflows!

Headquarters:

Chemspeed Technologies AG • Rheinstrasse 32 • 4302 Augst (BL) • Switzerland • Phone: +41 (0) 61 816 9500 • Fax: +41 (0) 61 816 9500
chemspeed@chemspeed.com

UK Office:

Chemspeed Technologies Ltd. • Caxton House, Northampton Science Park, Kings Park Road, Moulton Park • Northampton, NN3 6LG • United Kingdom
 Phone: +44 (1604) 654 200 • Fax: +44 (1604) 654 201 • uk@chemspeed.com

US Office:

Chemspeed Technologies Inc. • 113 North Center Drive • North Brunswick, NJ 08902 • USA • Phone: +1 (732) 329 1225 • Fax: +1 (732) 329 1226
us@chemspeed.com